



### Get Perfect Support - Remain Active!

The new *Invacare* **XLT** is the ideal active, rigid chair for day to day use. Designed to be perfectly adjusted to fit the body the **XLT** offers its users a wide range of various levels of support, enabling and even improving independence and mobility.



**Yes, you can.®**



### Top Configurability made Simple!

The *Invacare XLT* has been designed with the individual in mind. Measured, constructed and fitted to suit the precise requirements of each user, the **XLT** is the ideal chair to help regain, retain and improve both mobility and independence.

Complex chair-specific solutions and spare parts are a thing of the past with the new **XLT**. Benefit from the extensive and readily available range of Invacare options which can be fitted to the **XLT** without any hassle.

### Stability, Rigidity and Security for an Active Life

Made entirely from titanium the **XLT** gives its user stability and rigidity without compromising on weight. In combination with the wide range of options, features and configuration possibilities offered with the **XLT** the user will benefit from a chair which offers the maximum security and support needed when going about busy day to day life.

### Proven Quality

*Invacare*, together with *Rea*™ have made no compromises when it comes to the new **XLT**. Designed and built at the Rea factory in Sweden, users will benefit

from the strong tradition of both *Invacare* and *Rea* quality fused to create one ultimately unique chair: The *Invacare XLT*.

### Features and Options



#### Armrests

Wide range of armrests for best support.



#### Anti tipper device

The Anti-tipper device can easily be folded away by the foot. The device is also angled outwards, to give space for the feet for the attendant.



#### Standard wheel

with black Marathon Evolution Plus tyres.

**Features and Options**



**Front Caster Attachments**  
for stability and performance.



**Titanium Foottube**  
to improve comfort as well as lightness.

**Swing Away**  
for easy transfer and maneuverability as well as comfort.



**Tyres**  
① solid, ② pneumatic profiled or puncture proof, ③ pneumatic non marking, ④ low profile, ⑤ pneumatic profiled or puncture proof, ⑥ black marathon evolution.



**Frontwheels**  
① 75 x 24 hard (3"), ② 100 x 27 medium hard (4"), ③ 125 x 27 medium hard (5"), ④ 140 x 37 soft (5,5"), ⑤ 150 x 27 medium hard (6"), ⑥ 150 x 30 soft (6"), ⑦ 150 x 32 pneumatic (6"), ⑧ 180 x 40 pneumatic (7"), ⑨ 180 x 40 soft (7")

## Technical data



<b>Invacare XLT Active</b>	355 - 480 mm	400 / 450 / 500	400 - 520 mm	200 - 490 mm	380 - 505 mm	0 - 14°
<b>Invacare XLT Dynamic</b>	355 - 500 mm	400 / 450 / 500	400 - 500 mm	200 - 490 mm	355 - 485 mm	0 - 14°
<b>Invacare XLT Swing</b>	330 - 505 mm	400 / 450 / 500	390 - 520 mm	200 - 490 mm	280 - 510 mm	0 - 14°

in increments of 15 mm on 40 mm adj.



<b>Invacare XLT Active</b>	± 10°	SW + 200 mm	610 - 1050 mm	830 - 1030 mm	9.8 kg	135 kg
<b>Invacare XLT Dynamic</b>	± 10°	SW + 200 mm	610 - 1050 mm	730 - 880 mm	9.8 kg	135 kg
<b>Invacare XLT Swing</b>	± 10°	SW + 200 mm	610 - 1050 mm	820 - 1230 mm	11.3 kg	135 kg



<b>Invacare XLT Active</b>	6 kg
<b>Invacare XLT Dynamic</b>	6 kg
<b>Invacare XLT Swing</b>	6 kg

Crash-tested  
ISO 7176-19

The **Invacare XLT** complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: [www.invacare-rea.com](http://www.invacare-rea.com) under Rea wheelchairs/ Technical downloads.

## Frame colours



Blue



Red



Black



Grey



Blasted Titanium\*  
not for Swing



Orange\*

\* Available as special orders

## Upholstery



Black Nylon



Black Jemima

Invacare International Sàrl  
Route de Cité-Ouest 2 - 1196 Gland  
Switzerland  
Tel: +41 (0)22 354 60 10  
Fax: +41 (0)22 354 60 11  
[hqeuropa@invacare.com](mailto:hqeuropa@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

© 2011 Invacare International Sàrl  
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

XLT - EU - 10/2011



**Yes, you can.®**